

Tabletop Autoclave ITA-424



The ITA-424 front loading 'Benchtop' autoclave is pulse vacuum type and range is available in single sizes of 24 litres chamber capacity. These units require no plumbing, and only a standard single phase power supply. Build in 7 Sterilization cycles + 2 User definable cycles + 2 Test programs.

Technical Specification:

Model	ITA-424
Capacity (liters)	24 L
Chamber Dimensions (mm)	260x455 mm
Sterilisation Temperature (°C)	121°C (1.1 kg/cm²) / 134°C (2.1 kg/cm²)
External Dimensions (w x d x h) mm	560x665x400 mm
Net Weight (kg)	57 kg
Gross Weight (kg)	60 kg
Power (W)	2000 W
Electrical Requirements	AC 230V±10% 50HZ / 60HZ, 9A
Catalogue No.	580150502

Technical Data:

Model	ITA-424
Volume of Chamber	24 L
Air Expulsion	Pulse Vacuum
Sterilization Temp	121? (1.1 kg/cm²) / 134? (2.1 kg/cm²)
Dry Function	Adjustable 3 to 30 minutes
Water Reservoir	Supply tank : 4.5L / Discharge tank : 2.5L
Sterilization Programs	7 Sterilization cycles + 2 User definable cycles + 2 Test programs
Standard Accessory	Sterilization frame, sterilization network disk (375x220x20mm x 4 drive, Note: 3Kg / plate) SUS304 sterilization clip x1 branch network disk removed
Optional Accessory	Printer sets, SD Card
Power Consumption	2000 W
Net Weight	57 kg
Gross Weight	60 kg
Size of Chamber	260x455 mm
External Dimensions	560x665x400 mm
Electrical Requirement	AC 230V±10% 50HZ / 60HZ, 9A
Catalogue No.	580150502

Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States