

Benchtop Autoclave/Steam Sterilisers ITA-745



Benttop Autoclave/Steam Sterilisers ITA-745 provide the broadest range of steam sterilizers and has a solution available to maximize department productivity and optimize workflow. It can be used to Inactivated AIDS (HV), hepatitis B (HBV) mad cow virus and spores. The pressure door is auto locking device provide easy option for such a cost effective autoclave. Design with water supply and detection system that is used to clean disinfection with functions of water quality detection and debris filter. Mainly used in Dental clinic, medical, laboratory, Beauty Shop applications. Equipped with excellent features contain relief valve, pneumatic locking, and breakdown self-checking system with LCD display of pressure and temperature

Technical Specification:

Model	ITA-745
Capacity (liters)	45L
Chamber Dimensions (mm)	319x620 mm
Sterilisation Temperature (°C)	121°C, 134°C
Number of Shelves	2-layer movable tray
External Dimensions (w x d x h) mm	610x840x580 mm
Net Weight (kg)	100 kg
Power (W)	3500
Electrical Requirements	AC 220V/50Hz
Catalogue No.	580300502

Technical Data:

Model	ITA-745
Chamber	45L
Sterilization Temperature	121°C, 134°C
Tray	2-layer movable tray
Level	European Class B Standard
Safety	Safety valve, pneumatic lock, fault self-checking system
Test	Bowie & Dick steam penetration test, helix spiral test, Vacuum test, water quality testing
Drying program	Powerful vacuum drying, the remaining humidity of the instrument is <0.2%
Display	LCD
Applications	Inactivated AIDS (HV), hepatitis B (HBV) mad cow virus and spores
Interface	printer and USB
Power Consumption	3500 VA
External Dimension	610x840x580 mm
Internal Dimension	319x620 mm
Net Weight	100 kg
Electrical Requirements	AC 220V/50Hz
Catalogue Number	580300502

Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States